

Optical Tweezers: Measuring the Piconewton Forces



PHYS5381 Final Presentation

**Runze Shi; April 3th, 2026
Department of Physics
College of Science
Northeastern University**

Motivation


History, Mechanisms, Application.

- First demonstrated by Arthur Ashkin in 1986
- Essential in biophysics Research.
 - Elasticity of DNA
 - Stepping of molecular motors
- Induced force near the focal point due to intensity gradient.
- Restores for small displacement from the equilibrium position.



Material and Methods


The experimental setup



Material and Methods

The experimental setup


- Bead: trapped particles
 - Diameter of $4.16 \pm 0.21 \mu\text{m}$
 - Index of refraction $n = 1.50$
 - Deionized water $n = 1.33$



Material and Methods


Data Acquisition

- Measurements were taken at
 - 25, 50, 75, 120, 125, 130, 135, 145 $\mu\text{L}/\text{min}$.
- Python with the OpenCV, NumPy, and SciPy libraries.
- A region of interest (ROI): rows 240-270 was selected, as it trapped bead and optical trap center




Material and Methods

Data Acquisition



Results





Data analysis



Results

Data analysis


- Median Filter
- IQR Method
- Z-Score
- Combined



Results


Data analysis

- Savitzky-Golay filter was applied with
 - 46 frames (~1.5 s)
 - polynomial order 3
- Velocity vs. time computed from taking the derivative of Position vs. time
- Negative velocity minima is selected.



Results

Savitzky-Golay filter



Results

Data analysis




- Viscous drag force a trapped bead was calculated using Stokes' law
 - η : the dynamic viscosity of water at room temperature
 - r : bead radius
 - v : is the measured free-bead velocity at each flow rate
- The trap stiffness k was determined by fitting
 - a linear model to the measured force-displacement pairs across all flow rates:

$$F_{\text{drag}} = 6\pi\eta r v$$

$$F = kx$$

Displacement & Force vs. Flow Rate

- (ODR) fit was applied to the force displacement data
- 1 sigma error band shown in red around ODR.
- Displacement is $\sim 2 \mu\text{m}$ and drag force is $\sim 0.065 \text{ pN}$ at the max flow rate $\sim 140 \mu\text{L/min}$



Conclusion

- Successfully calibrated a single-beam optical trap using viscous drag from controlled flow rates
- Trap stiffness measured with the $R^2 = 0.713$; Confirmed Hookean behavior over 0 - 3 μm displacement range
- Maximum sustainable trapping force $\sim 0.065\text{pN}$, consistent with the laser power and NA used.
- All data and analysis pipeline were all uploaded in the GitHub repository <https://github.com/Lemon-Elf/PHYS5318.git>; branch optical_tweezers